

## THE NU ATTOM HIGH RESOLUTION ICP-MS: SEMI-QUANTITATIVE ANALYSIS FOR RAPID SAMPLE SCREENING

### INTRODUCTION

Fast, multi-elemental analysis for preliminary evaluation of sample concentrations is a must in routine laboratory work. This kind of screening has been proven to be a useful tool to help choose the best analytical technique or as a “ballpark” figure to determine approximate concentration figures for further analytical investigation (i.e. determination of required dilution factors, concentration calibration standard range, or for spiking for ID-TIMS analysis).

Quadrupole ICP-MS instruments are typically used for this application. However, high resolution ICP-MS (HR-ICP-MS) instruments may also have a role to play due to their superior detection limits, better sensitivity and ability to resolve isobaric interferences.

In this note we describe the use of a Nu Attom HR-ICP-MS in a semi-quantitative screening role to determine its suitability for this application.



### Instrumentation

The Nu Attom is a double-focusing, high-resolution magnetic sector mass spectrometer. The instrument is entirely purpose designed and built to provide the best performance and reliability coupled with flexibility and ease-of-use for precise and accurate elemental and isotope ratio analysis. A unique detector system gives the Nu Attom a large dynamic range, and its electrostatic scanning capability has the widest range in its class (40%). Furthermore, the continuously variable high resolution means that sufficient resolution for isobaric separation can be selected with minimum compromise in sensitivity.

### Experiment

The first step in this experiment was to generate a mass response curve using a multi-elemental standard (SPEX CertiPrep Group, USA) diluted to 1 ppb using 2% HNO<sub>3</sub> and spiked with an internal standard. Blanks, standards and samples were all spiked in this way with 1 ppb Rh for internal normalization processes. Internal standard normalization and blank correction were applied to the mass response curve. The standard used to determine the mass response curve contained only a selection of elements / isotopes across the mass range and not necessarily every target element / isotope. They are listed below:

Li(7)  
Al(27)  
Sc(45)  
Mn(55)  
Y(89)  
In(115)  
Cs(133)  
Tl(205)  
U(238)

Once the mass response curve had been established, an unknown sample was analysed and the concentrations were calculated directly from this curve.

A sample of the USGS standard material BCR-2 (Basalt, Columbia River, USA) was treated as the “unknown” sample and measured using the semi-quantitative mode of the Attom HR-ICP-MS.

The BCR-2 sample was taken from a stock solution and was diluted by a factor of 10 with 2% HNO<sub>3</sub> for analysis. A blank consisting of 2% HNO<sub>3</sub> was also prepared. The 1 ppb Rh spike was added to the blank and sample solution. No mathematical corrections for isobaric overlap (e.g. Sn on Cd) or spectral interferences (e.g. BaO on Eu) were made.

### Discussion

The results are reported in Table 1 and consist of 10 repeat analysis of the standard BCR-2. Data are shown along with certified, recommended or previously reported concentrations wherever available.

In addition, Table 1 shows the accuracy values (in the far right column) represented by the percentage difference between the measured value and the certified (or recommended / previously reported) value. It can be seen that all the results are within  $\pm 20\%$  of the recommended values for this particular reference material.

### Conclusions

The Nu Attom is a high resolution ICP-MS that is ideal for the most demanding ICP-MS requirements where sensitivity, precision, and speed of analysis are paramount. This note highlights that the instrument is capable of being used in a rapid screening mode for the semi-quantitative analysis of multi-element samples.

BCR-2 (n=10)	ppm					
	Average	Stdev	RSD%	Certified	Stdev	Difference%
Be(9)	0.3	0.0	2.9	N/A		
Sc(45)	32.5	0.4	1.1	33.0	2	-1.5
V(51)	406.0	3.1	0.8	416.0	14	-2.4
Mn(55)	1749.8	11.6	0.7	1520.0	60	15.1
Co(59)	32.8	0.8	2.3	37.0	3	-11.3
Sr(88)	361.7	2.7	0.7	340.0	3	6.4
Y(89)	34.6	0.2	0.7	37.0	2	-6.6
Zr(90)	182.9	1.3	0.7	184.0	1	-0.6
Nb(93)	10.9	0.1	0.5	12.6	0.4	-13.4
Mo(95)	252.4	1.5	0.6	250.0	20	1.0
Ag(107)	0.2	0.0	11.2	N/A		
Ag(109)	0.1	0.0	5.8	N/A		
Cd(111)	0.1	0.0	10.5	N/A		
In(115)	0.2	0.0	9.1	N/A		
Sn(119)	2.1	0.1	3.7	N/A		
Sn(120)	2.1	0.1	2.6	N/A		
Sb(121)	0.1	0.0	5.2	N/A		
Cs(133)	1.0	0.0	2.1	1.1	0.1	-9.0
Ba(137)	615.3	5.3	0.9	677.0	2	-9.1
La(139)	24.0	0.3	1.1	24.9	0.2	-3.6
Ce(140)	53.3	0.5	0.9	52.9	0.2	0.8
Pr(141)	7.3	0.1	1.4	6.7	0.1	9.2
Nd(143)	31.9	0.5	1.5	28.7	0.1	11.0
Nd(145)	32.3	0.4	1.1	28.7	0.1	12.5
Nd(146)	32.4	0.3	1.0	28.7	0.1	12.8
Sm(147)	7.4	0.1	1.6	6.6	0.02	12.4
Sm(149)	7.5	0.2	2.0	6.6	0.02	13.9
Eu(153)	2.4	0.0	1.8	2.0	0.01	20.4
Gd(157)	7.6	0.2	2.2	6.8	0.03	13.0
Tb(159)	1.2	0.0	1.7	1.1	0.03	15.3
Dy(163)	7.7	0.1	1.4	6.4	0.05	20.0
Ho(165)	1.6	0.0	1.8	1.3	0.03	21.2
Er(166)	4.4	0.1	1.9	3.7	0.01	20.8
Tm(169)	0.6	0.0	3.1	0.5	0.04	18.3
Yb(172)	4.2	0.1	2.8	3.4	0.02	24.0
Lu(175)	0.6	0.0	3.0	0.5	0.009	17.3
Hf(178)	5.7	0.1	1.9	4.9	0.1	15.5
Ta(181)	0.8	0.0	3.2	0.7	0.02	4.4
W(182)	0.5	0.0	3.4	N/A		
Re(185)	0.0	0.0	27.8	0.01	0.003	-17.4
Hg(202)	0.1	0.0	13.4	N/A		
Tl(205)	0.3	0.0	4.6	N/A		
Pb(208)	9.4	0.2	2.5	11.0	1	-14.3
Bi(209)	0.0	0.0	15.3	N/A		
Th(232)	5.8	0.1	2.0	5.7	0.5	1.2
U(238)	1.8	0.0	1.9	1.7	0.19	9.2

Table 1 : The data displayed is from 10 repeats of BCR-2 (Basalt, Columbia River, USA). Also displayed are certified, recommended or previously reported values where these are available. The far right hand column shows the differences between the measured and certified values.