



ASTRUM HIGH RESOLUTION GD-MS

Performance Specifications

Resolution:	Continuously variable resolution
Average Background:	<0.5 counts per second (cps)
Operating Analyser Pressure:	<2 x 10 ⁻⁷ mbar
Matrix Current:	>5x10 ⁻¹⁰ amps at RP > 4000 (Using total Cu signal)
Signal Stability:	<3% RSD over 10 minutes (Using Cu @ RP=4000) <5% RSD over 30 minutes (Using Cu @ RP=4000)
Peak Side Stability:	<+/-30ppm over 30 minutes (Using mass 63)
Acquisition Mass Stability:	<100ppm from Li to U
Ion Counting Efficiency:	>70% relative to Faraday
Abundance Sensitivity:	mass 62 (baseline) / Cu63 <200ppb at RP > 4000
Internal Reproducibility:	Major < 3%RSD Minor < 5%RSD Trace < 10%RSD
External Reproducibility:	Major < 5%RSD Minor < 10%RSD Trace < 20%RSD
Gas Background:	C < 1ppm (gas and matrix RSF values set to 1) N < 1ppm (gas and matrix RSF values set to 1) O < 1ppm (gas and matrix RSF values set to 1)



ASTRUM

HIGH RESOLUTION GD-MS

Nu Instruments introduce the next generation of glow discharge mass spectrometers. The ion source has been purpose designed in collaboration with users to provide the best in performance, reliability and flexibility.

The double-focusing, high-resolution magnetic sector mass spectrometer is based on the successful Nu AttoM HR-ICP-MS for accurate and precise elemental analysis.

The no-compromise design of the Nu Astrum GD-MS provides the best in solids analysis.

Resolution

Continuously variable resolution. Both Source and Collector slits are continuously variable and computer controlled

Average Background

<0.5 counts per second (cps), Average baseline signal at masses where no ion beam signal exists, for example 220.5 amu.

Operating Analyser Pressure

< 2×10^{-7} mbar, under full gas load. The analyser section of the instrument (ESA, Flight Tube, and Collector) is all metal sealed and pumped using ion pumps

Matrix Current

> 5×10^{-10} amps at RP > 4000, using total Cu signal of a Cu sample

Signal Stability

Using a Cu acquisition routine with parameters of RP > 4000, 200 points per peak window, 2 sweeps per isotope, and 1 cycle, the peak area of masses 63 and 65 using a Cu sample will be:

<3% RSD over 10 minutes (5 repeats)

<5% RSD over 30 minutes (15 repeats)

Peak Side Stability

At RP > 4000, the peak side stability will be shown to be <+/-30ppm over 30 minutes using a Cu sample and monitoring mass 63

Acquisition Mass Stability

Using the data acquisition program, the position of peaks will be shown to be within 100ppm in mass (demonstrated using peaks from C up to Ta)

Ion Counting Efficiency

The electron multiplier will have an efficiency of >70% relative to the Faraday amplifier when the multiplier supply voltage is set within the plateau

Abundance Sensitivity

Using a Cu matrix and a resolution of > 4000, the abundance sensitivity measured at mass 62 will be <200ppb relative to mass 63

Internal Reproducibility

Using a homogeneous reference material, measurements will be made on one major (approx 10 – 50%), one minor (approx 0.1%), and one trace element (approx 1ppm). From five repeat measurements, the following reproducibility will be demonstrated:

Major < 3%RSD

Minor < 5%RSD

Trace < 10%RSD

External Reproducibility

The sample used in the previous test will be withdrawn, dismantled then re-inserted and a further five analyses performed. The difference between the certified concentrations and the average concentrations measured will be:

Major < 5%RSD

Minor < 10%RSD

Trace < 20%RSD

Gas Background

Using a pure Ta sample, and by measuring ion beam ratio relative to Ta, the levels of gas background will be shown to be less than:

C < 1ppm

N < 1ppm

O < 1ppm

Assuming that the RSF values for the gases and the matrix element(s) are set to 1.000. Measurements will be made after a suitable pre-sputter time using a cryo-cooled cell